

# Lentivirus HIV-1 p24 Instructions

**Cat:EHY0002**

## Content

	CAT	Volume
1 CP (Coated Plate)	EHY0002CP	96 well
2 S (Standard)	EHY0002S1~S7,S0	8 vial
3 DA-H (Detect Antibody-HRP)	EHY0002DA-H	6 ml/bottle
4 SD (Sample Diluent)	ESD01	15 ml/bottle
5 VL (Virus Lysis)	EVL01	12 ml/bottle
6 TS (TMB Substrate)	ETS01	12 ml/bottle
7 SS (Stop Solution)	ESS01	12 ml/bottle
8 WB (Wash Buffer 10×)	EWB01	50 ml/bottle
9 SF (Sealer Film)	ESF01	6 pieces

**NOTE:** After the kit is opened, the stabilization period of each content is 30 days, so please use the kit within 30 days after opening.

## Sample Dilution

Samples of the crude virus or filtered need a 20-fold dilution into Sample Diluent. A suggested 20-fold is 10  $\mu$ l of sample + 190  $\mu$ l of Sample Diluent.

For chromatographic samples, at least 100-fold dilution is required. A suggested 100-fold is 5  $\mu$ l of sample + 495  $\mu$ l of Sample Diluent.

## REAGENT PREPARATION

### ■ Washing Buffer (1×) Preparation

Pour entire contents (50 ml) of the **Washing Buffer Concentrate** (10×) into a clean 500 ml graduated cylinder. Bring to final volume of 500 ml with glass-distilled or deionized water. Transfer to a clean wash bottle and store at 2 to 25°C.

### ■ Standard Curve Preparation:

S1 to S7 and S0 is ready to use.

## ASSAY PROCEDURE

**Bring all reagents and samples to room temperature before use.**

- 1 Prepare all reagents and working standards as directed in the previous sections.
- 2 Remove excess **CP** (Coated Plate) strips from the plate frame, return them to the foil pouch and reseal.
- 3 Add 50  $\mu$ l of **VL** (Virus Lysis) to each well.
- 4 Add 10  $\mu$ l of **Standard or sample** per well. Ensure reagent addition is uninterrupted and completed within 15 minutes.
- 5 Add 50  $\mu$ l of **DA-H** (Detect Antibody-HRP) to each well.
- 6 Cover with an **SF** (Sealer Film). Incubate at room temperature (18 to 25°C) for 30 mins on a microplate **shaker** set at 500 rpm.
- 7 Aspirate each well and **wash**, repeating the process four times. Wash by filling each well with **WB** (Washing Buffer 300  $\mu$ l). Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining **WB** (Washing Buffer) by aspirating or decanting. Invert the plate and **blot** it against clean paper towels.
- 8 Add 100  $\mu$ l of **TS** (TMB Substrate) to each well. Incubate for 5-30 minutes at room temperature.
- 9 Add 100  $\mu$ l of **SS** (Stop Solution) to each well.
- 10 Determine the optical density within 30 minutes, using microplate **reader** set to 450 nm corrected with 570 nm or 630 nm.

The minimum detectable dose (MDD) of HIV-1 p24 is typically less than 0.12 ng/ml.

The MDD was determined by adding two standard deviations to the mean optical density value of ten zero standard replicates and calculating the corresponding concentration.


## PRECISION

**Intra-assay Precision (Precision within an assay)** Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.

**Inter-assay Precision (Precision between assays)**

Sample Number	Intra-assay Precision			Inter-assay Precision		
	S1	S2	S3	S1	S2	S3
	22	22	22	6	6	6
Average (ng/ml)	62.5	333.8	928.4	62.6	324.3	978.5
Standard deviation	3.1	15.9	58.3	1.8	6.0	39.4
Coefficient of variation (%)	4.9	4.8	6.3	3.9	1.9	4.0

## TYPICAL DATA



ng/ml	O.D.	Average	Corrected
0.00	0.007	0.006	
	2	8	0.0070
	0.022	0.021	
4.12	0	2	0.0216
	0.045	0.046	0.0146
12.35	3	1	0.0457
	0.108	0.106	0.0387
37.04	5	6	0.1076
	0.220	0.218	0.1006

## SENSITIVITY

## RECOVERY

The spike recovery was evaluated by spiking 3 levels of HIV-1 p24. The recovery ranged from 95% to 106% with an overall mean recovery of 101%.

## LINEARITY

To assess the linearity of the assay, five samples were spiked with high concentration of HIV-1 p24. The linearity ranged from 97% to 103% with an overall mean recovery of 100%.